



**Marble Medical**  
Your ideas. Our capabilities. Your product.

## Clear Protective Face Shields

Affordable protection for healthcare<sup>1</sup>, sanitation, industrial, and food service professionals.



**ALL MATERIALS SOURCED AND FABRICATED IN THE U.S.A.**

**DESCRIPTION:** Face shield designed for protection from liquids or dust.  
ONE SIZE FITS ALL  
10mil Polystyrene Shield – 10.5” tall x 11.75” wide  
Anti-Fog Coating – Nonsterile – Not Shatter Proof – Latex Free  
Open Cell Foam Forehead Cushion  
One-Wrap Hook & Loop Adjustable Strap

**PRODUCTION CAPACITY:** 30,000 Pieces per week

**DIRECTIONS:** Peel liner from foam & apply foam to top of shield  
Feed hook & loop strap through each slot – run across front of shield  
Place Face Shield on head  
Adjust the Strap to fit foam comfortably above the eyebrows.

**FOR INQUIRIES OR ORDERS:** Call us at (850) 385-4441  
Email us at [info@marblemed.com](mailto:info@marblemed.com)

**PRICING: \$2.10 EACH - MINIMUM ORDER QUANTITY – 160 PIECES**

---

<sup>1</sup> This face shield is intended for single use by a single user for medical purposes. The product has not been cleared or approved by the FDA. The FDA has authorized the use of this product by healthcare providers as personal protective equipment (PPE) under an Emergency Use Authorization (EUA) for the duration of the EUA. <https://www.fda.gov/medical-devices/emergency-situations-medical-devices/emergency-use-authorizations#covid19ppe>,