



Marble Medical
Your ideas. Our capabilities. Your product.

Clear Flexible Protective Face Shield Sheets

Inexpensive contamination protection for general consumer use.



Clothing provided by Knot Responsible®

ALL MATERIALS SOURCED AND FABRICATED IN THE U.S.A.

DESCRIPTION: Face shield sheet designed for protection from liquids or dust particles.

ONE SIZE FITS ALL

1mil Polypropylene Shield – 24” tall x 18” wide

Nonsterile – Not Shatter Proof – Latex Free

1/2” Tape Strips at Top and Bottom of Sheet

PRODUCTION CAPACITY: 20,000 Pieces per week

DIRECTIONS: Peel one paper liner from tape strip & adhere to a cap or visor above the bill

Peel the remaining liner and adhere tape strip to clothing above the chest

FOR INQUIRIES OR ORDERS: Call us at (850) 385-4441

Email us at info@marblemed.com

PRICING: \$.35 EACH - MINIMUM ORDER QUANTITY – 2,000 SHEETS

VOLUME PRICING AVAILABLE UPON REQUEST