



Marble Medical
Your ideas. Our capabilities. Your product.

One Piece Clear Protective Face Shields

Affordable protection for healthcare¹, sanitation, industrial, and food service professionals.



ALL MATERIALS SOURCED AND FABRICATED IN THE U.S.A.

DESCRIPTION: Face shield designed for protection from liquids or dust.
ONE SIZE FITS MOST (20" to 25 1/8" circumference)
Dimension of shield 9" tall x 13" wide
Open Cell Foam Forehead Cushion
Anti-Fog Coating – Nonsterile – Not Shatter Proof – Latex Free
Protective film prevents scratching while in transit
200 units per carton (24" x 24"x 14")

PRODUCTION CAPACITY:10,000 Pieces per week

DIRECTIONS: Remove protective film from front of shield
Place Face Shield on head
Adjust the foam forehead cushion to fit comfortably above the eyebrows.

FOR INQUIRIES OR ORDERS: Call us at (850) 385-4441
Email us at info@marblemed.com

PRICING: \$1.90 EACH – VOLUME PRICING AVAILABLE
MINIMUM ORDER QUANTITY – 200 PIECES

¹ This face shield is intended for single use by a single user for medical purposes. The product has not been cleared or approved by the FDA. The FDA has authorized the use of this product by healthcare providers as personal protective equipment (PPE) under an Emergency Use Authorization (EUA) for the duration of the EUA. <https://www.fda.gov/medical-devices/emergency-situations-medical-devices/emergency-use-authorizations#covid19ppe>